



Clinical Multitest Chemistry Analyzer BAS-120TS

TOUCH SCREEN Multitest Chemistry Analyzer

**THIS BIOCHEMISTRY ANALYZER
CAN TEST SAMPLES IN TWO WAYS,
BY SUCTION PUMP AND BY
DISPOSABLE CUVETTE**



*The BAS-120TS is open to being used
with all brands of biochemistry reagents.*

Clinical Biochemistry Touch Screen Analyzer with Incubator Temperature Control and Internal Printer BAS-120TS is a TOUCH SCREEN Microprocessor Biochemical Analyzer, which measures over 50 biological indicators.

Features

- **New Version**
- **Powerful Software**
- **AD Auto Zero**
- **7.0" Touch Screen**
- **Lamp Standby Function**
- **Multi-Language**
- **Auto Temperature Control**
- **Built-in Printer**
- **24V Universal Power Supply**
- **USB Port for external printer, online operation, keyboard, software update (optional)**
- **Very easy to use, user-friendly interface, requiring minimal training or supervision to get started.**
- **50 Biochemistry tests**

Accessories

- **1 power cable**
- **1 pump cassette**
- **1 waste tube**
- **1 aspiration tube**
- **2 rolls of printer paper**
- **1 touch pen**
- **1 lamp house**
- **5 clotting cuvettes**
- **1 clotting cuvette adapter**

Test Menu

T. BILI	Total Bilirubin	Ca	Calcium
D. BILI	Direct Bilirubin	CL	Chlorine
TP	Total Protein	P	Phosphorus
ALB	Albumin	C	Carbon
Lactic Acid	Lactic Acid	A-HBDH	alpha-hydroxy butyrate dehydrogenase
LP (a)	Lipoprotein A	AMY	Amylase
CHE	Cholinesterase	ADA	Adenosine Deaminase
Fe	Iron	C3	Complement Component 3
IgM	Immunoglobulin M	APOB-IB	Apolipoprotein-B
R-GT/GGT	Reflocheck Glucose/Gamma Glutamyl Transferase	APOB	Apolipoprotein-B
ALP/AKP	Alkaline Phosphatase	CK-MB	Creatine Kinase-MB
UREA	Urea	ASO	Antibody Streptolysin O
CREA	Creatinine	CO2	Carbon Dioxide
ALT/GPT	Alanine Aminotransferase/Glutamate Pyruvate Transaminase	LDL-C	
AST/GOT	Aspartate Aminotransferase/Serum glutamic-oxaloacetic transaminase	HS-CRP	High Sensitivity C-Reactive Protein
LA		C4	Complement Component 4
Cu	Copper	CHE	Acetylcholinesterase
IgA	Immunoglobulin A	Mg	Magnesium
CHO	Chinese Hamster Ovary	AFU	Alpha-L-Fucosidase
TG	Thyroglobulin	FMN	Flavin Mono Nucleotide
CK	Creatinine Kinase	HDL-C	High Density Lipoprotein-C
LDH	Lactate Dihydrogenase	PA	
UA	Urinalysis	IgA	Immunoglobulin A
GLU	Glucose		
TBA	Total Bile Acids		
Zn	Zinc		

and drug toxicity