## <u>hybe</u>x

## Media Storage Bottles



hybex<sup>™</sup> reusable media bottles are manufactured from high quality, borosilicate 3.3 glass. With their superior chemical resistance, these bottles are ideal for storage of reagents, culture media, biological fluids and a variety of other aqueous and nonaqueous solutions.

A low coefficient of thermal expansion (3.3), provides temperature resistance that is far superior to that of standard lab glass (or soda lime glass). Hybex glass can withstand temperatures from -70 to 500°C, permitting repeated use in freezers, microwaves, autoclaves and dry heat sterilization chambers.

All hybex bottles have permanent, white enamel graduations for volume measurement and a white identification box for marking or coding.

All hybex media bottles include a drip-free sealing ring as well as a standard blue polypropylene (GL45) sealing cap. Bottles are also available with green, purple, white and yellow caps, or in an assortment pack containing two of each color.

- **Autoclavable caps (GL45) & drip rings**
- Permanent enamel graduation markings
- Large area for bottle marking/identification
- Autoclavable, reuseable & breakage resistant
- Five cap colors

## Technical Data

Material: hybex<sup>™</sup> borosilicate glass (3.3)

Cap Thread Size: GL45\*

Temp. Range: -70°C to 500°C Sterilization: Autoclave or dry heat Cap Material:

Polypropylene (140°C) Cap Colors: Blue, green, purple

yellow, white or assorted

\*50ml bottles are only supplied with blue GL32 sealing caps

## Ordering Information

B3000-RIN2

B3000-50 Media storage bottle, 50ml with standard (GL32) blue cap, 10/pk. Media storage bottle, 100ml with standard (GL45) blue cap, 10/pk. Media storage bottle, 250ml with standard (GL45) blue cap, 10/pk. B3000-100\* B3000-250\* B3000-500\* Media storage bottle, 500ml with standard (GL45) blue cap, 10/pk. Media storage bottle, 1L with standard (GL45) blue cap, 10/pk. Media storage bottle, 2L with standard (GL45) blue cap, 5/pk. B3000-1000\* B3000-2000\* B3000-SP\* Media storage bottle, *Starter Pack* with standard (GL45) blue cap (Starter Pack of ten bottles includes 2x100ml, 3x250ml, 3x500ml & 2x1000ml)

B3000-CAP\* High temperature red caps, up to 180°C (GL45), 10/pk B3000-HTC Replacement cap, standard blue (GL45), 10/pk. B3000-RIN Replacement sealing ring, (GL45),10/pk. Replacement cap, standard blue (GL32), 10/pk. B3000-CAP2

Replacement sealing ring, (GL32), 10/pk. \*Please add (-G) for green caps, (-P) for purple caps, (-Y) for yellow caps, (-W) for white caps or (-AS) for assorted caps (2 of ea. color)



Benchmark Scientific Inc. **PO Box 709** Edison, NJ 08818 Phone: 908-769-5555

Fax: 908-222-1864

Web: www.BenchmarkScientific.com Email: Info@BenchmarkScientific.com