


Benchmark 


INCU-SHAKER *Series*

Shaking Incubators

- + *Exceptional temperature uniformity*
- + *Miniature and cooling options available*
- + *Exchange flask clamps in one second with the new **MAGiC CLAMP™** magnetic system*



INCU-SHAKER Mini




INCU-SHAKER 10L

WITH
MAGiC CLAMP






Incu-Shaker 10L/LR
(H2010/H2012)



Incu-Shaker Mini
(H1001-M)



H1010-P-MP



H1010-P-SP



H1000-MR-1550

Technical Data:

Temperature Range (Mini):	Ambient +5 to 70°C (Increments of 0.1°C)
Temperature Range (10L):	Ambient +5 to 60°C (Increments of 0.1°C)
Temperature Range (10LR):	Ambient -15 to 60°C (Increments of 0.1°C)
Temperature Uniformity:	+/- 0.25 %
Temperature Accuracy:	+/- 0.1 %
Shaking Speed Range:	30 to 300 rpm (Increments of 1 rpm)
Shaking Orbit:	19 mm (3/4")
Platform Dimensions:	
Incu-Shaker Mini	9.5 x 11.5 in. / 23.5 x 29.2 cm
Incu-Shaker 10L/LR	17.8 x 17.8 / 45 x 45 cm
Exterior Dimensions:	
Incu-Shaker Mini	11 x 18.8 x 13 in. / 28 x 48 x 33 cm
Incu-Shaker 10L	20 x 26 x 19 in. / 51 x 65 x 47 cm
Incu-Shaker 10LR	20 x 31 x 22 in. / 51 X 79 X 56 cm.
Weight:	
Incu-Shaker Mini	25 lbs. / 11.5 kg
Incu-Shaker 10L	120 lbs / 55 kg
Incu-Shaker 10LR	170 lbs / 78 kg
Electrical:	120 or 230V, 50-60Hz

Ordering Information

H1001-M* Incu-Shaker Mini, includes standard rubber mat platform, 115V
 H2010* Incu-Shaker 10L, includes standard rubber mat platform, 115V
 H2012* Incu-Shaker 10LR, includes standard rubber mat platform, 115V

**To order an item in 230V, add (-E) to item number.*

Optional Platforms and Shelves

H1000-MR MAGic Clamp universal platform, for Incu-Shaker Mini (9.5 x 11.5")
 H1010-MR MAGic Clamp universal platform, for Incu-Shaker 10L/LR (18 x 18")
 H1000-P-MP Microplate platform, 4 position for Incu-Shaker Mini
 H1010-P-MP Microplate platform, 6 position for Incu-Shaker 10L/LR
 H1000-P-SP Spring rack platform, holds uncommon vessels for Incu-Shaker Mini
 H1010-P-SP Spring rack platform, holds uncommon vessels for Incu-Shaker 10L/LR
 H1010-SH Optional Shelf, for Petri Dishes, microplates, trays, etc. (14 x 5.5") for Incu-Shaker 10L/LR (max. 3)

Flask & Tube Rack Accessories for MAGic Clamp Platforms

H1000-MR-50 Flask Clamp, 50 ml (H1000-MR: max. 20 / H1010-MR: max. 49)
 H1000-MR-125 Flask Clamp, 125 ml (H1000-MR: max. 12 / H1010-MR: max. 36)
 H1000-MR-250 Flask Clamp, 250 ml (H1000-MR: max. 8 / H1010-MR: max. 25)
 H1000-MR-500 Flask Clamp, 500 ml (H1000-MR: max. 5 / H1010-MR: max. 16)
 H1000-MR-1000 Flask Clamp, 1000 ml (H1010-MR: max. 9)
 H1000-MR-2000 Flask Clamp, 2000 ml (H1010-MR: max. 5)
 H1000-MR-MP Microplate Clamp, (H1000-MR: max. 4 / H1010-MR: max. 8)
 H1000-MR-1550 Tube Rack, 20x50 & 30x15 ml tubes (H1000-MR: max.1 / H1010-MR: max. 4)
 H1000-MR-T15 Tube Rack, 32x15 tubes (H1000-MR: max. 2 / H1010-MR: max. 8)
 H1000-MR-T50 Tube Rack, 15x50 tubes (H1000-MR: max. 2 / H1010-MR: max. 8)
 H1000-MR-T600 Bottle Rack, 4x500/600ml conical bottles (H1000-MR: max. 1 / H1010-MR: max. 2)
 H1000-MR-TSR MAGic Clamp™ tilted holder for 15 ml and 50 ml tube racks
 H1000-MR-CMB Adjustable flask clamp for MAGic clamp platform, accepts 125 ml (max. 36), 250 ml (max. 25) or 500 ml (max. 16)

Additional accessories available. Please visit www.BenchmarkScientific.com for more details.

Maximum capacity of Erlenmeyer flasks

Flask Size:	Incu-Shaker Mini	Incu-Shaker 10L/LR
50ml	20	49
125ml	12	36
250ml	8	25
500ml	5	16
1L	2	9
2L	0	5



See it on the web!