





LABOMED, INC.

www.labomed.com
spectro@labomed.com



Bullet Holder

Filters

Polarizing Attachment

LB-1010 Comparison Microscope

LB-1010 Comparison Microscope can observe two objects with a pair of eyepiece at the same time. Using field cutting, jointing and overlapping methods, two (or more two) objects can be compared together. LB-1010 has clear image, high resolution, and can identify tiny differences between objects accurately. It is basically used for forensic science and police departments. LB-1010 has excellent optical performance and complete comparison function, it is suitable for various comparison demands, so it has a very wide range of applications.