

Standard Features

- LCD display
- Quick response
- Accumulation, counting and percentage modes
- 13 different units of measure
- Detachable enclosed weighing chamber (in capacities ≤1 kg)
- Battery operation for up to 40 hours of continuous use (without backlight)
- RS-232 interface for connecting to printer or PC
- Plastic enclosure

Specifications

Weighing Units

Gram, Ounce, Pound, Carat, Pennweight, Troy Ounce, Grain (See technical manual for complete listing)

Operating Modes

Weighing, Percent Weighing, Parts Counting, Accumulation

Power

AC adapter (12 V/500 mA); internal rechargeable battery (6 V/1.2 Ah)

Stabilization Time

2 seconds

Operating Temperature

32° F to 104° F (0° C to 40° C)

Calibration Weight

EQ-300: 300 g, EQ-600: 600 g, EQ-1500: 1,500 g
EQ-3000: 3,000 g, EQ-6000: 6,000 g

Est. Weight

6.6 lb (3 kg)

Warranty

One-year limited warranty

Return Policy

Rice Lake Standard Policy

Battery Life

Up to 40 hours of continuous use (without backlight)

Approval

CE Marked



Model	Capacity	Readability	Repeatability	Linearity	Weighing Chamber	Pan Size	Part #	List Price
EQ-300	300 g	0.01g	0.01g	0.02 g	6.7 x 5.7 x 3.1in (170 x 145 x 80 mm)	Ø 4.7 in (Ø 120 mm)	165338	Consult
EQ-600	600 g	0.02 g	0.02 g	0.04 g	6.7 x 5.7 x 3.1in (170 x 145 x 80 mm)	Ø 4.7 in (Ø 120 mm)	165339	Consult
EQ-1500	1,500 g	0.05 g	0.05 g	0.1g	N/A	4.8 x 5.6 in (122 x 142 mm)	165340	Consult
EQ-3000	3,000 g	0.1g	0.1g	0.2 g	N/A	4.8 x 5.6 in (122 x 142 mm)	165341	Consult
EQ-6000	6,000 g	0.2 g	0.2 g	0.4 g	N/A	4.8 x 5.6 in (122 x 142 mm)	165342	Consult

Options/Accessories

Part#	Description	Price
41279	Printer cable, DB9M to 25-Pin M	Consult
103142	SP742 dot matrix printer, tear bar	Consult
119919	AbleCable RS-232/USB cable (requires 130789)	Consult
130789	Slim line null modem adapter	Consult
131135	Aluminum case	Consult
165183	AC adapter, 110-240 VAC, interchangeable plugs	Consult



PN 165183 Multi-range AC adapter and multi-plug

