

Standard Features

- External weight calibration
- Large, 16.5 mm-high backlit liquid crystal display (LCD)
- Mono-metal tuning fork sensor (MMTS) provides quick response and stability
- 20-step bar graph display
- One part in 620,000 resolution
- RS-232 interface for connecting to printer or PC
- Plastic housing on models <1,000 g
- Plastic top and die cast bottom on models >1,000 g

OPTIONS

- Rechargeable battery
- Under-hook weighing

Specifications

Weighing Units

Carat, Pennyweight, Troy Ounce, Decimal Ounce, Decimal Pound (See technical manual for complete listing)

Operating Modes

Weight, Count, Percentage, Limit Function (Checkweighing)

Power

120 VAC (adapter included)

Battery Option

Nickel Cadmium battery, 12-hour run time, 24-hour charge time

RS-232 Baud Rate

1200/2400/2800/9600; Data bits - 8, Connector - DIN 5-Pin

Calibration Weight

TP-220: 200 g
TP-820: 800 g
TP-3200: 3,000 g
TP-8200: 8,000 g
TP-12K: 12,000 g

Operating Temperature

41° F to 95° F (5° C to 35° C)

Est. Weight

6.2 lb (2.8 kg)

Warranty

Five-year limited warranty

Return Policy

Rice Lake Standard Policy

Model	Capacity	Readability	Repeatability	Linearity	Pan Size	Weighing Chamber	Dimensions	Part #	List Price
TP-220	220 g	0.001 g	1 division	1 division	Ø 4.6 in (116.8 mm)	170mm W x 140 mm D x 95 mm H	9.2 x 7.1 x 6.6 in (including windshield) (235 x 182 x 170 mm)	107244	Consult
TP-820	820 g	0.01 g	1 division	1 division	6.6 in x 5.5 in (167.6 x 139.7 mm)	N/A	9.2 x 7.1 x 2.9 in (235 x 182 x 75 mm)	107245	Consult
TP-3200	3,200 g	0.01 g	1 division	1 division	7 in x 6.2 in (177.8 x 157.5 mm)	N/A	10.1 x 7.5 x 3.4 in (265 x 192 x 87 mm)	107246	Consult
TP-8200	8,200 g	0.1 g	1 division	1 division	7 in x 6.2 in (177.8 x 157.5 mm)	N/A	10.1 x 7.5 x 3.4 in (265 x 192 x 87 mm)	107247	Consult
TP-12K	12,000 g	0.1 g	1 division	1 division	7 in x 6.2 in (177.8 x 157.5 mm)	N/A	10.1 x 7.5 x 3.4 in (265 x 192 x 87 mm)	107248	Consult

Options/Accessories

Part #	Description	Price
107249	Rechargeable battery ≤TP-820	Consult
107237	Rechargeable battery ≥TP-3200	Consult
107250	In-use dust cover ≥TP-3200	Consult
107240	In-use dust cover ≤TP-820	Consult
109236	Underweigh hook ≤TP-820	Consult
109292	Underweigh hook ≥TP-3200	Consult
42566	115 VAC adapter/12 VDC	Consult
42168	230 VAC adapter	Consult
107251	Additional operating manual	Consult
75351	WinWedge® Tal software	Consult
115303	5-pin to 9-pin DB female RS-232 cable	Consult

